

## Declaration on Widex DAILY

To whom it may concern

- All WIDEX DAILY hearing aids in the various form factors and in all price point 50, 100, 110, 220, 330 and 440 are electro-acoustical identical and share the same hardware platform with the similar form factors for the generic WIDEX DREAM hearing aid as well as the specific price points 440, 330, 220 and 110
- “MPO, acoustical gain and technical settings such as EMC immunity are not changed and are already measured and stated in the current Delta reports for the generic WIDEX DREAM form factor covered by the following reports:

	Covering the <u>generic</u> version and the price points:	Delta reports for: Acoustical Tests IEC 60118-0,1,2,6	Delta reports for: EMC IEC 60118-13 ANSI C63.19	Delta reports for: Australian Acoustic requirements, OHS DoSO
DREAM D-PA	440, 330, 220, 110	DREAM D-PA S: PR036_12  DREAM D-PA M: PR037_12	DREAM D-PA S: EMC036_12  DREAM D-PA M: EMC037_12  DREAM D-PA S: Ansi036_12  DREAM D-PA M: Ansi037_12	DREAM D-PA S: PR217_12  DREAM D-PA M: PR218_12
DREAM D-FS	440, 330, 220, 110	DREAM D-FS S BTE: PR137_12  DREAM D-FS M BTE: PR138_12  DREAM D-FS P BTE: PR135_12  DREAM D-FS HP BTE: PR136_12	DREAM D-FS S: EMC137_12  DREAM D-FS M: EMC138_12  DREAM D-FS P: EMC135_12  DREAM D-FS HP: EMC136_12	DREAM D-FS S: NR137_12  DREAM D-FS M: NR138_12  DREAM D-FS P : NR135_12A  DREAM D-FS HP: NR136_12A

			DREAM D-FS S: Ansi137_12  DREAM D-FS M: Ansi138_12  DREAM D-FS P: ANSI135_12  DREAM D-FS HP: Ansi136_12	
DREAM D-9	440, 330, 220, 110	DREAM D-9: PR214_11	DREAM D-9: EMC214_11  DREAM D-9: Ansi214_11	DREAM D-9: PR215_12
DREAM D-FA	440, 330, 220, 110	DREAM D-FA: PR078_13	DREAM D-FA: EMC078_13  DREAM D-FA: Ansi078_13	DREAM D-FA: PR079_13
DREAM D-FA P	440, 330, 220, 110	DREAM D-FA P: PR195_14	DREAM D-FA P: EMC195_14  DREAM D-FA P: ANSI195_14	DREAM D-FA P: PR197_14
DREAM D-XP	440, 330, 220, 110	DREAM D-XP: PR168_12	DREAM D-XP: EMC168_12  DREAM D-XP: Ansi168_12	DREAM D-XP: PR225_12
DREAM D-CIC	440, 330, 220, 110	DREAM D-CIC: PR076_12  DREAM D-CIC TR: PR100_12	DREAM D-CIC: EMC076_12  DREAM D-CIC TR: EMC100_12  DREAM D-CIC: Ansi076_12  DREAM D-CIC TR: Ansi100_12	DREAM D-CIC: PR219_12  DREAM D-CIC TR: PR220_12

DREAM D-CIC-M	440, 330, 220, 110	Reports N/A. Ref. to. "Declaration on Widex DREAM CIC-M and DREAM CIC-M TR	DREAM D-CIC-M: EMC100_14 DREAM D-CIC-M TR: EMC101_14 DREAM D-CIC-M: ANSI100_14 DREAM D-CIC-M TR: ANSI101_14	Reports N/A. Ref. to. "Declaration on Widex DREAM CIC-M and DREAM CIC-M TR
---------------	--------------------	---	--	---



02-02-2015 Ole Gelsbo  
Date In-market Product Mgr, Global Marketing



02-02-2015 Hans-Otto Bindeballe  
Date Certification Mgr. Regulatory affairs